## High resolution spectroscopy without cross-dispersion

### Guillem Anglada-Escude

University of Hertfordshire & Queen Mary University of London guillem.anglada@gmail.com

#### Ulrike Lemke

University of Goettingen ulrike.lemke@phys.uni-goettingen.de

## Outline of the talk

The stabilized echelle spectrometer


Cross-dispersion using wavelength sensitive pixels

Feasibility using current technology


Trade-offs over classic design

Key challenges


## The stabilized echelle spectrometer




## The stabilized echelle spectrometer




## The stabilized echelle spectrometer









**HARPS spectrometer**European Southern Observatory – La Silla






Carnegie's PFS spectrometer Credits: J Crane. OCIW



**Order curvature** 




Line tilt and optical aberrations

Optimal 2D extraction nightmarish Modeling the spectra on the detector might be possible but... no success, requires perfectly known:

- > PSF maps
- > Wavelength maps
- > True stellar template

Still get correlated flux measurements if the spectra is to be analyzed.


# Cross dispersion using wavelength sensitive pixels



We only need to assign a photon to a diffraction order

Grove size is the only key parameter for the echelle grating

resolution by the post-disperser needed



Ricardson Gratings


ARCONS: A 2024 PIXEL OPTICAL THROUGH NEAR-IR CRYOGENIC IMAGING SPECTROPHOTOMETER

#### Mazin et al. 2013 PASP

#### **MKIDS**

Microwave Kinetic Inductance Detectors

Status: 10 000 element array being produced (Ben Mazin Lab, UCSB)



The typical spectral resolution for the MKIDs used in ARCONS is R~10 at 4000 Å declining linearly with in-

The achieved spectral resolution is drastically below the theoretical limit for a 100 mK operating temperature of R=100 at 4000 Å (Mazin et al. 2012). We believe the

Measure phase of the microwaves coming from a resonant cavity

ARCONS: A 2024 PIXEL OPTICAL THROUGH NEAR-IR CRYOGENIC IMAGING SPECTROPHOTOMETER

#### Mazin et al. 2013 PASP

#### **MKIDS**

Microwave Kinetic Inductance Detectors

Status: 10 000 element array produced and tested on-sky (Ben Mazin, UCSB)


The typical spectral resolution for the MKIDs used in ARCONS is R~10 at 4000 Å declining linearly with in-

The achieved spectral resolution is drastically below the theoretical limit for a 100 mK operating temperature of R=100 at 4000 Å (Mazin et al. 2012). We believe the

#### Theoretically limited to order sorting of 40 A


3800-6800 A **>** 75 orders


HARPS 72 orders x 4096 pixels



Use IFU technology?

Image credits : Ulrike Lemke





#### VIRUS: a massively-replicated IFU spectrograph for HET

Gary J. Hill<sup>a\*</sup>, Phillip J. MacQueen<sup>a</sup>, Carlos Tejada<sup>b</sup>, Francisco J. Cobos<sup>b</sup>, Povilas Palunas<sup>a</sup>, Karl Gebhardt<sup>c</sup>, and Niv Drorv<sup>c</sup>

<sup>a</sup>McDonald Observatory, University of Texas at Austin, 1 University Station, C1402, Austin, TX 78712 -0259, USA


bInstituto de Astronomía de la Universidad Nacional Autónoma de México, Apdo. Postal 70-264, 04510 México

Department of Astronomy, University of Texas at Austin, 1 University Station, C1402, Austin, TX 78712 -0259 USA

#### ABSTRACT

We investigate the role of industrial replication in the construction of the next generation of spectrographs for large telescopes. In this paradigm, a simple base spectrograph unit is replicated to provide multiplex advantage, while the engineering costs are amortized over many copies. We argue that this is a cost-effective approach when compared to traditional spectrograph design, where each instrument is essentially a one-off prototype with heavy expenditure on engineering effort. As an example of massive replication, we present the design of, and the science drivers for, the Visible IFU Replicable Ultra-cheap Spectrograph (VIRUS). This instrument is made up of 132 individually small and

#### VIRUS HET →500\$ per pixel




http://hetdex.org/hetdex/virus.php

# Objective Lens #1 F = 800 mm WI-separation of 7 Angstrom FMOS-IFU Objective Lens #2 Low dispersion element Different orders (7 nm apart ?)


#### **Minispectrometers**


Optical layout in mini-spectrometer (RC series)



Hammamatsu catalog

#### **Optical**





Satisfy all design requirements!

Range 350-1150 nm


Resolution for a 50micron fiber : 16 A  $\rightarrow$ 

500 echelle orders

http://www.stellarnet.us/p



#### **Near infrared**



Same, but with InGaAs 1D array 10 000 EUR/unit

Very simple optical layout





Image credits : Ulrike Lemke



## Trade-offs over classic design

#### The promise

Very simple stabilized optics
No optical aberrations, order tilting or spectral line tilts
Detector decoupled from stabilized body
Detectors upgradeable at any time
No stability requirements for detectors
Very fast read-out, minimal read-out noise
Off-the-shelf technology
Minimal stray-light
Ultra-broad wavelength range

#### The challenges

2000 EUR x 4000 → 8 MEUR Fabrication of ultrastable 1D fiber coupling Performance of mini-spectrometers unknown

#### Unsolved issues

Large echelle grating → Slicing? Reformating? Input stabilitzation : scrambling

Funding proposals being submitted, collaborations are VERY WELCOME

## Trade-offs over classic design

#### The promise

Very simple stabilized optics

No optical aberrations, order tilting or spectral line tilts

Detector decoupled from stabilized body

Detectors upgradeable at any time

No stability requirements for detectors

Very fast read-out, minimal read-out noise

Off-the-shelf technology

Minimal stray-light

Ultra-broad wavelength range

Ask Ulrike Lemke, ulrike.lemke@phys.uni-goettingen.de Fiber input feed and detailed optic designs

The challenges

2000 EUR x 4000 → 8 MEUR

Fabrication of ultrastable 1D fiber coupling Performance of mini-spectrometers unknown

Guillem Anglada-Escude guillem.anglada@gmail.com Overall instrument concept and echelle design

#### Unsolved issues

Large echelle grating → Slicing? Reformating? Input stabilitzation : scrambling


Funding proposals being submitted, collaborations are VERY WELCOME

## Summary

- The post-multiplexed spectrometer can be built today
- Many advantages over classic echelle design
- All technology available. Very few key elements need testing

   → mini-spectrometers, 1D fiber array couplings
- Upgradeable detectors
- Project on-hold, funding pending

Candidate technology for a optical+nlR ultrastabilized ELT-HIRES spectrometer for the E-ELT?



**Upgradeable**